

Quantitative Rapid Testing for Mycotoxins in Cereals



10min lateral flow assay

dipstick rapid testing in 10min

Test 5 mycotoxins with same preparation

afatoxin, T2 toxin, vomitoxin, fumunisin, zearalenone as well as ochratoxin


Works for for creals, feed

such as maize, wheat and rice

* Quantitative result with NBReader® Fast 2

NBReader Fast 2 can be used to determine the result

Sample Preparation



Result determination



Insert the ID card to reset the standard curve.

Insert the test strip into the NBReader® Fast 2.

Check the result.

* The test result can be determined visually by eye if NBReader is not available



Scan to read PDF copy

RING BIOTECHNOLOGY CO LTD





Ringbio

Quantitative Rapid Testing for Mycotoxins in Cereals

Mycotoxins Rapid Test Kit Technical Specifications

Product #	Product Name	Detection range (ug/kg, ppb)
Q401	IQuanti Aflatoxin B1 Rapid Test Kit	0 - 20 or 0 - 30
Q601	IQuanti Zearalenone Rapid Test Kit	0 - 150
Q602	IQuanti Vomitoxin Rapid Test Kit	0 - 4000
Q603	IQuanti Fumonisin Rapid Test Kit	0 - 5000
Q605	IQuanti Ochratoxin Rapid Test Kit	0 - 40
Q604	IQuanti T2-toxin Rapid Test Kit	0 - 800
M1003	NBReader ® Fast 2 system	/

We also have qualitative rapid tests please contact us for more information.

Other feed test kits



Product #	Product Name
F1001	Aflatoxin ELISA Kit
F1002	Zearalenone ELISA Kit
F1003	Deoxynivalenol ELISA Kit
E6006	Melamine ELISA Kit
E6002	Aflatoxin M1 ELISA Kit

You can always visit www.ringbio.com for our latest product update

Not all test kits are registered in all countries, consult your sales before order

For veterinary use only, for research use only, for outside of China only.



Ring Biotechnology Co Ltd

Add: Building 3, Zhongtongtai TechnoPark, BDA, Beijing 100176, China

Tel: +86-10-56267496 / +86-10-26242030

Email: info@ringbio.com

Web: www.ringbio.com